

Cost effective solutions for office-based procedures

OfficePACK Includes:

- Disposable Mayo Tray: 1
- Mayo Tray Drape: 1
- Plastic Cotton Tip: 2
- Gauze 100x100mm: 15
- Surgical Lubricant 3g: 1
- Povidone Iodine Solution w/
3 Plastic Cotton Tip: 1
- Alcohol Prep Pad: 1
- Spinal Needle 22 GA: 1
- Sanitary Paper:
- Cysto Irrigation Set In flow w/
Luer Lock Fitting: 1
- Outflow Tubing w/ Luer Lock
Fitting: 1
- Gloves: 2
- OfficeSPEC Disposable
Side-Opening Speculum w/
Light: 1
- Endoscopic Seal 1.2mm w/
Luer Lock Fitting: 1
- Item Tray: 1
- Under Buttocks Drape w/
Graded Collection Bag: 1



Office The Premier Procedure Kit for Office Hysteroscopy **PACK**[™]

OBP Medical is proud to offer the in-office physician the OfficePACK Procedure Kit for Hysteroscopy. With Several exciting features never before seen in an all-in-one pack including: a custom mayo drape with built in instrument pockets and trash container, OfficeSPEC disposable side-opening speculum with light, 1.2mm Endoscopic double sealing seal, inflow and outflow tubing with luer-lock adaptors, under buttocks drape with graded drain bag, and a drawstring for easy

post procedure disposal. All Features exclusive to the OBP Medical OfficePACK Procedure Kit for Hysteroscopy. See below for a full list of the pack's contents.

Product Features:

- Includes OfficeSPEC disposable side-opening speculum w/ built in light source and 1.2mm Endoscopic Seal.
- Inflow and outflow tubing fitted with luer-lock adaptors compatible with any feamile luer-lock fitting.
- Custom mayo drape with built-in instrument pockets and trash container.
- Mayo tray drape doubles as trash bag with ties for easy post procedure clean up.
- All items come packed sterile on a disposable mayo tray ready for use.

ORDERING INFORMATION

- Item # C040100
- Sold in boxes of 5
- All items come sterile and are latex free.
- Call OBP Medical for additional product information or visit www.obpmedical.com.

OBP
MEDICAL

OBP Medical Inc.
360 Merrimack St Building 9
Lawrence, MA 01843

1-888-300-2946

OBP Medical is very proud to be an ISO13485 certified company.

For further information please visit www.obpmedical.com or call 1-888-300-2946