



M E D I C A L

Products for Office-Based Procedures

# Endoscopic SEAL

Leak Free Endoscopy

The OBP Medical Self Sealing Endoscopic Seal for Urology is available in two sizes:



0.6mm (Yellow Band) accepts instruments 2-5fr



1.6mm (Green Band) accepts instruments 4-10fr

*The OBP Medical Endoscopic Seal for Urology is used to prevent fluid leakage or air aspersation through the working channel of a cystoscope or ureteroscope during introduction and use of an instrument. All OBP Medical Endoscopic Seals feature a funnel guided top for easy insertion of instruments, self sealing technology for leak free use, and a luer-lock adaptor compatable with any female luer-lock fitting.*

- Funnel guided top for easy insertion of instrument
- Double seal technology for maximum leak control
- Self sealing for leak free use with multiple instruments
- Luer-lock adapter for use with any female luer-lock fitting



OBP MEDICAL ENDOSCOPIC SEALS FOR UROLOGY ARE SOLD IN PACKS OF 10 FOR \$105.00/ PER PACK

1.6MM SEAL 10 PACK ITEM # C010005    0.6MM SEAL 10 PACK ITEM # C010006

FOR MORE INFORMATION PLEASE VISIT [WWW.OBPMEDICAL.COM](http://WWW.OBPMEDICAL.COM) OR CALL 1-888-300-2946

**OBP is proud to provide the medical community with high quality medical devices and innovative solutions to address existing and emerging needs within the medical market. OBP is ISO 13485 certified to ensure production of only the finest in medical garments and equipment.**